

## Preparation of cells for FACS analysis: Farnham lab protocol; modified 02/26/08

### Prepare cells (this can be done up to two weeks in advance)

1. Trypsinize adherent cells using standard protocol (non-adherent cells can be counted directly).
2. Determine cell number. Use  $10^5$ - $10^6$  cells/sample.
3. Add appropriate volume of cells to a conical tube and spin them down at 1000 rpm for 3 min. Aspirate off the media.
4. Vortex pellet at low speed, add 0.5 ml of cold PBS, vortex again for 2-3 seconds. ***It is very important to achieve a single cell suspension.*** If cell clumping is a problem, you can use calcium- and magnesium-free PBS.
5. Resuspend the pellet in 5 ml of cold PBS. Centrifuge cells for 6 min at 1000 rpm.
6. Aspirate off the PBS.
7. Add 0.5 ml cold PBS and pipette up and down to achieve a single cell suspension.
8. Prepare a tube with 4.5 ml ice-cold 100% ethanol. Begin vortexing the tube with ethanol as you slowly add 0.5 ml of cells in PBS, dropwise.
9. Incubate on ice or at 4°C for at least 30 minutes. At this point, samples can be stored in the refrigerator for up to two weeks.
10. Make an appointment with the Flow Cytometry Facility!

### Stain cells (this is done immediately before FACS analysis)

1. Prepare propidium iodide (PI)/Triton staining solution with RNase A.  
Recipe: 10 ml of 0.1% (v/v) Triton X-100 (Sigma) of PBS (Triton is a viscous liquid), 2 mg DNase-free RNase A (Sigma), 200 microliters of 1 mg/ml PI (Molecular Probes). Prepare fresh!
2. Warm tubes to 37°C for 5-10 minutes (otherwise ice crystals can tear apart your cells).
3. Centrifuge the ethanol-suspended cell for 5 min at 200 g (=1000 rpm). Decant ethanol thoroughly.
4. Suspend cell pellet in 5 ml of PBS, wait about 1 min, centrifuge 5 min at 200 g (=1000 rpm) for 5 min. Discard the supernatant.
5. Suspend cell pellet in 1 ml PI/Triton X-100 staining solution with RNase A. Keep in the dark at RT for 30 min or at 37°C for 15 min.
6. Immediately before analysis, vortex each sample and filter using a 5 ml polystyrene round-bottom tube with cell-strainer cap (Becton Dickinson, cat# 352235). Pour cell suspension onto the cap and centrifuge the tubes for a short period of time at 1000 rpm.
7. You are ready to do the FACS analysis.